

# Alzheimer's Disease

## Clinical utility of A $\beta$ 42/40 ratio



What is the clinical utility of plasma A $\beta$ 42/40 ratio in patients being evaluated for Alzheimer's disease?



### Background


One of the characteristics of Alzheimer's disease is the accumulation of beta-amyloid (A $\beta$ ) plaques in the brain, which have traditionally been assessed by positron emission tomography (PET) imaging or measurement of cerebrospinal fluid. While the role of blood-based biomarkers in the assessment of A $\beta$  plaques is still being explored, blood-based biomarkers, including the ratio of peptides A $\beta$ 42/40 in plasma, have been found to correspond with A $\beta$ -PET status and may guide or complement PET.



### Methods and Results

The study population included 250 participants who provided plasma samples and underwent an amyloid PET scan. A $\beta$ 42/40 ratio was measured by liquid chromatography-tandem mass spectrometry (LC-MS/MS).

#### Diagnostic performance of plasma A $\beta$ 42/40 ratio for A $\beta$ -PET status<sup>1</sup>



Patients with A $\beta$ 42/40 ratio  $\geq 0.170$  have a low likelihood of A $\beta$ -PET positivity (Alzheimer's disease pathology) and may be able to avoid PET imaging, thus allowing doctors to focus on other causes of cognitive impairment.

1. Weber DM, Taylor SW, Lagier RJ, et al. Clinical utility of plasma A $\beta$ 42/40 ratio by LC-MS/MS in Alzheimer's disease assessment. *Front Neurol.* 2024;15:1364658. doi:10.3389/fneur.2024.1364658

For more information about Quest Diagnostics publications, please visit the [Clinical Education Center](#).