

Liver Fibrosis, Fibrosis-4 (FIB-4) Index Panel

Test Code: 30555, 30710

Specimen Requirements: Refrigerated serum (1 mL in serum separator tube [SST®], 0.5 mL minimum); **and** room-temperature whole blood (full EDTA lavender-top tube, 1 mL minimum; or full EDTA lavender Microtainer®, 0.5 mL minimum)

CPT Codes*: 84450; 84460; 85049

CLINICAL USE

• Assess likelihood of advanced liver fibrosis in patients with nonalcoholic fatty liver disease (NAFLD), hepatitis B, or hepatitis C

CLINICAL BACKGROUND

NAFLD, hepatitis B, and hepatitis C are liver diseases that can result in abnormal development of scar tissue (ie, fibrosis). NAFLD is considered a manifestation of metabolic syndrome within the liver and is associated with type 2 diabetes and cardiovascular disease,^{1,2} whereas hepatitis B and C are caused by infection with the hepatitis B (HBV) and C (HCV) viruses, respectively.^{3,4} Many patients with these liver diseases have a mild clinical presentation. However, chronic disease progression can lead to the development of fibrosis, cirrhosis, and associated complications such as hepatocellular carcinoma.¹⁻⁴

Management of NAFLD, hepatitis B, and hepatitis C depends, in part, on liver disease severity that is characterized by the presence of inflammation and/or advanced fibrosis.⁵⁻⁷ Tests for advanced fibrosis are of particular interest owing to the strong association of fibrosis stage with mortality in NAFLD⁸ and management implications for hepatitis B and hepatitis C.^{6,7} A biopsy is the gold standard for assessing liver disease severity, including the presence of fibrosis, but the procedure is highly invasive and involves safety risks.³⁻⁵ Biopsies are also susceptible to sampling error (ie, a representative tissue sample is not collected).⁵ Thus, noninvasive tests that assess for advanced fibrosis have been developed.

One such noninvasive test, the FIB-4 index, uses a combination of routine blood tests to indicate whether a patient has a high or low probability of advanced fibrosis.

The FIB-4 index formula⁹ yields a single score by combining patient age with measurements of 3 biomarkers: aspartate aminotransferase (AST), alanine aminotransferase (ALT), and platelet count.

The FIB-4 index has clinical utility for management of NAFLD, hepatitis B, and hepatitis C, and is included in respective clinical practice guidance:

- For NAFLD, the American Association for the Study of Liver Diseases (AASLD) recommends the FIB-4 index as an option for noninvasive testing to identify patients with high likelihood of advanced fibrosis; patients with an elevated FIB-4 index are recommended to consider further evaluation with a biopsy.⁵ The FIB-4 index has a high negative predictive value (≥90%) in patients with NAFLD, which makes this test particularly useful for ruling out advanced fibrosis and thereby avoiding unnecessary biopsies.^{10,11} In addition, this test is useful in a primary care setting to determine whether a patient with NAFLD requires referral to a specialist for further evaluation.
- For hepatitis B, AASLD guidance includes the FIB-4 index as an alternative to biopsy for assessing disease severity; this assessment helps guide treatment decisions for patients with chronic HBV infection, including identification of appropriate candidates for antiviral therapy and determination of therapy duration.⁶
- For hepatitis C, a joint panel from AASLD and the Infectious Diseases Society of America recommends the FIB-4 index as an option for noninvasive testing to evaluate individuals with HCV infection for advanced fibrosis.⁷ This guidance considers an elevated FIB-4 index to be evidence of cirrhosis, which, in turn, determines treatment strategy selection and indicates additional surveillance for associated complications.⁷

Quest Diagnostics offers the Liver Fibrosis, Fibrosis-4 (FIB-4) Index Panel (test code 30555), which includes tests for AST, ALT, and platelet count, plus calculation of the FIB-4 index. Quest also offers the Liver Fibrosis, Hepatic Function Panel with Fibrosis-4 (FIB-4) Index (test code 30710), which includes the FIB-4 index in conjunction with additional tests for liver function: total protein, albumin, globulin, albumin/ globulin ratio, bilirubin (total, direct, and indirect), and alkaline phosphatase.



Table. Interpretation Information for the FIB-4 Index, by Testing Indication¹²

Indication for testing	Compatible with the absence of advanced fibrosis ^a	Compatible with the presence of advanced fibrosis ^a	Indeterminate result
NAFLD	<1.30	>2.67	1.30-2.67
Hepatitis B	<1.00	>2.65	1.00-2.65
Hepatitis C	<1.45	>3.25	1.45-3.25

NAFLD, nonalcoholic fatty liver disease.

^aAdvanced fibrosis is defined as stage 3 or stage 4.

INDIVIDUALS SUITABLE FOR TESTING

- Individuals with NAFLD who are undergoing risk stratification and will potentially be considered for a liver biopsy
- Individuals with chronic HBV infection who are being considered for, or receiving, antiviral therapy
- Individuals with HCV infection

METHOD

- AST and ALT measured using spectrophotometry
- Platelet count measured using electronic cell sizing and counting/cytometry/microscopy
- The FIB-4 index calculated with a formula using AST, ALT, platelet count, and patient age

INTERPRETIVE INFORMATION

The Table contains information useful in interpreting FIB-4 index values for NAFLD, hepatitis B, and hepatitis C. For all testing indications, an elevated FIB-4 index is consistent with the presence of advanced fibrosis, whereas a low FIB-4 index is consistent with the absence of advanced fibrosis. Patient characteristics and clinical features should guide interpretation. Consider additional assessment for patients with an indeterminant result.

References

- 1. Marjot T, Moolla A, Cobbold JF, et al. Nonalcoholic fatty liver disease in adults: current concepts in etiology, outcomes, and management. *Endocr Rev.* 2020;41(1):66-117. doi:10.1210/ endrev/bnz009
- Stefan N, Häring HU, Cusi K. Non-alcoholic fatty liver disease: causes, diagnosis, cardiometabolic consequences, and treatment strategies. *Lancet Diabetes Endocrinol*. 2019;7(4):313-324. doi:10.1016/s2213-8587(18)30154-2

- Tang LSY, Covert E, Wilson E, et al. Chronic hepatitis B infection: a review. JAMA. 2018;319(17):1802-1813. doi:10.1001/ jama.2018.3795
- 4. Spearman CW, Dusheiko GM, Hellard M, et al. Hepatitis C. *Lancet*. 2019;394(10207):1451-1466. doi:10.1016/s0140-6736(19)32320-7
- Chalasani N, Younossi Z, Lavine JE, et al. The diagnosis and management of nonalcoholic fatty liver disease: practice guidance from the American Association for the Study of Liver Diseases. *Hepatology*. 2018;67(1):328-357. doi:10.1002/ hep.29367
- Terrault NA, Lok ASF, McMahon BJ, et al. Update on prevention, diagnosis, and treatment of chronic hepatitis B: AASLD 2018 hepatitis B guidance. *Hepatology*. 2018;67(4):1560-1599. doi:10.1002/hep.29800
- Ghany MG, Morgan TR. Hepatitis C guidance 2019 update: American Association for the Study of Liver Diseases-Infectious Diseases Society of America recommendations for testing, managing, and treating hepatitis C virus infection. *Hepatology*. 2020;71(2):686-721. doi:10.1002/hep.31060
- Ekstedt M, Hagström H, Nasr P, et al. Fibrosis stage is the strongest predictor for disease-specific mortality in NAFLD after up to 33 years of follow-up. *Hepatology*. 2015;61(5):1547-1554. doi:10.1002/hep.27368
- Sterling RK. Fibrosis-4 (FIB-4) index for liver fibrosis. MDCalc. Accessed February 26, 2021. https://www.mdcalc.com/fibrosis-4fib-4-index-liver-fibrosis
- Bril F, McPhaul MJ, Caulfield MP, et al. Performance of plasma biomarkers and diagnostic panels for nonalcoholic steatohepatitis and advanced fibrosis in patients with type 2 diabetes. *Diabetes Care*. 2020;43(2):290-297. doi:10.2337/dc19-1071
- Shah AG, Lydecker A, Murray K, et al. Comparison of noninvasive markers of fibrosis in patients with nonalcoholic fatty liver disease. *Clin Gastroenterol Hepatol*. 2009;7(10):1104-1112. doi:10.1016/j.cgh.2009.05.033
- 12. Kumar R, Teo EK, How CH, et al. A practical clinical approach to liver fibrosis. *Singapore Med J*. 2018;59(12):628-633. doi:10.11622/smedj.2018145
- *The CPT codes provided are based on AMA guidelines and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.

QuestDiagnostics.com

Quest, Quest Diagnostics, any associated logos, and all associated Quest Diagnostics registered or unregistered trademarks are the property of Quest Diagnostics. All third-party marks—[®] and TM—are the property of their respective owners. © 2021 Quest Diagnostics Incorporated. All rights reserved. TS9819 02/2021