

# ANCA Screen with MPO and PR3, with Reflex to ANCA Titer

**Test Code:** 70159(X)

**Specimen Requirements:** 2 mL room temperature serum; 0.8 mL minimum

**CPT Codes\*:** 86021 (x3)

## CLINICAL USE

- Differential diagnosis of systemic vasculitis
- Monitoring patients with systemic vasculitis

## CLINICAL BACKGROUND

Systemic vasculitis is characterized by inflammation of and damage to blood vessels. Subsequent disruption of the blood supply leads to tissue and organ damage. Vasculitis may be a primary or secondary manifestation of disease and can be caused by certain infections, malignancy, rheumatic disease, medications, and a wide variety of autoimmune disorders. Differential diagnosis is required to determine the appropriate therapy.

Differential diagnosis can be aided by testing for specific antineutrophil cytoplasmic antibodies (ANCA), which have been associated with several of the autoimmune systemic vasculitis disorders. These include granulomatosis with polyangiitis (GPA, formerly Wegener's), microscopic polyangiitis (MPA), and eosinophilic granulomatosis with polyangiitis (EGPA, formerly Churg-Strauss).<sup>1</sup> Each disorder is associated with predominance of a specific ANCA type.<sup>2</sup> The

ANCA types are revealed by fluorescent patterns obtained in the indirect immunofluorescence ANCA screen. For example, the cytoplasmic pattern (C-ANCA) is very common in GPA, but not MPA or EGPA. The perinuclear pattern (P-ANCA), on the other hand, is rare in GPA, common in MPA, and moderately common in EGPA cases. The atypical P-ANCA pattern is rare in all 3 of these; it is usually associated with nonvasculitic conditions such as inflammatory bowel disease.<sup>3</sup> The sensitivity and specificity of these markers for the various disorders are summarized in the **Table**.<sup>4-6</sup>

The diagnostic accuracy of the ANCA screen can be improved by combining it with immunoassays specific for myeloperoxidase (MPO) and proteinase-3 (PR3) antibodies. An international consensus group recommends this approach.<sup>2</sup> Though the C-ANCA pattern typically reflects specificity to PR3, there is not 100% concordance between C-ANCA and PR3 antibody as C-ANCA has multiple targets. Similarly, the P-ANCA pattern predominantly reflects MPO specificity.

ANCA can also be used to help guide patient management. A persistently positive ANCA is associated with relapse if treatment is stopped;<sup>7</sup> thus continued treatment should be considered as long as the ANCA is positive. ANCA titers may start high, decline during treatment, and increase again in relapse. ANCA titers do not closely reflect disease activity, though.<sup>7</sup> Thus when titers increase, the patient should be monitored more closely but not necessarily receive increased treatment.<sup>7</sup> Additional measures of disease activity should be used before modifying treatment.<sup>7,8</sup>

**Table. Diagnostic Sensitivity (Specificity) (%) of ANCA, PR3 Antibody, and MPO Antibody for Various Diseases**

Marker	Vasculitides				Other Autoimmune Disease			
	GPA <sup>a,6</sup>	MPA <sup>a,6</sup>	EGPA <sup>4</sup>	PAN <sup>a,6</sup>	GCA <sup>a,6</sup>	UC <sup>5</sup>	SLE <sup>a,6</sup>	RA <sup>a,6</sup>
ANCA	85 (93)	68 (87)	31	15	3	55	11	9
C-ANCA	81 (100)	3 (93)	5	0	1	ND	0	0
C-ANCA+/PR3+	69 (100)	0	1	0	2	ND	0	0
P-ANCA	4 (94)	65 (94)	21	15	2	55	11	9
P-ANCA+/MPO+	2 (99)	48 (100)	20	0	1	15	2	0

GPA, granulomatosis with polyangiitis; MPA, microscopic polyangiitis; EGPA, eosinophilic granulomatosis with polyangiitis; PAN, polyarteritis nodosa; GCA, giant cell arteritis; UC, ulcerative colitis; SLE, systemic lupus erythematosus; RA, rheumatoid arthritis; and ND, not determined.

<sup>a</sup> Sensitivity and specificity based on patients tested for ANCA in a rheumatology clinic.

## INDIVIDUALS SUITABLE FOR TESTING

- Individuals with symptoms of systemic vasculitis (eg, unexplained systemic illness with multiple organ involvement, unexplained ischemia)

## METHOD

### ANCA Screen

- Cell-based, indirect immunofluorescence using fixed human neutrophils
- Results reported: negative or positive
- Positive results reflexed to titer of the relevant pattern(s) (eg, C-ANCA, P-ANCA, atypical P-ANCA) at additional charge (CPT code 86021)

### Myeloperoxidase and Proteinase-3 Antibodies

- Semiquantitative multiplex immunoassay
- Analytical sensitivity: 0.2 AI for each antibody
- Analytical specificity: no known cross reactivity with antibodies associated with infection and immune disorders

Panel components can be ordered separately: ANCA screen (70171), MPO antibody (8796), PR3 antibody (34151).

## REFERENCE RANGES

ANCA Screen	Negative (titer <1:20)
Myeloperoxidase Antibody	<1.0 AI (no antibody detected)
Proteinase-3 Antibody	<1.0 AI (no antibody detected)

## INTERPRETIVE INFORMATION

A positive ANCA screen supports a diagnosis of autoimmune-related systemic vasculitis in a symptomatic patient (**Table**). Positive results are also seen in inflammatory bowel disease (ulcerative colitis) and occasionally in other autoimmune diseases (systemic lupus erythematosus, rheumatoid arthritis, autoimmune hepatitis). Exposure to certain drugs (eg, propylthiouracil, hydralazine, methimazole) and infectious agents (eg, hepatitis C virus) can result in secondary vasculitis and an ANCA-positive screen result.<sup>8,9</sup>

A negative ANCA, MPO antibody, and/or PR3 antibody result does not rule out systemic vasculitis.

An increase in serum ANCA levels suggests the need for closer clinical monitoring.<sup>7</sup> Patients who are persistently ANCA negative after treatment may be at less risk for a relapse.<sup>10</sup>

Owing to limitations in sensitivity and specificity, ANCA, MPO antibody, and PR3 antibody test results should be interpreted carefully in light of clinical and other laboratory data.

## References

- Jennette JC, Falk RJ, Bacon PA, et al. 2012 revised International Chapel Hill Consensus Conference Nomenclature of Vasculitides. *Arthritis Rheum.* 2013;65:1-11.
- Savigne J, Dimech W, Fritzler M, et al. Addendum to the International Consensus Statement on testing and reporting of antineutrophil cytoplasmic antibodies. Quality control guidelines, comments, and recommendations for testing in other autoimmune diseases. *Am J Clin Pathol.* 2003;120:312-318.
- Papp M, Altorjay I, Lakos G, et al. Evaluation of the combined application of ethanol-fixed and formaldehyde-fixed neutrophil substrates for identifying atypical perinuclear antineutrophil cytoplasmic antibodies in inflammatory bowel disease. *Clin Vaccine Immunol.* 2009;16:464-470.
- Comarmond C, Pagnoux C, Khellaf M, et al. Eosinophilic granulomatosis with polyangiitis (Churg-Strauss): clinical characteristics and long-term followup of the 383 patients enrolled in the French Vasculitis Study Group cohort. *Arthritis Rheum.* 2013;65:270-281.
- Reese GE, Constantinides VA, Simillis C, et al. Diagnostic precision of anti-Saccharomyces cerevisiae antibodies and perinuclear antineutrophil cytoplasmic antibodies in inflammatory bowel disease. *Am J Gastroenterol.* 2006;101:2410-2422.
- Schönermarck U, Lamprecht P, Csernok E, et al. Prevalence and spectrum of rheumatic diseases associated with proteinase 3-antineutrophil cytoplasmic antibodies (ANCA) and myeloperoxidase-ANCA. *Rheumatology (Oxford).* 2001;40:178-184.
- Ntatsaki E, Carruthers D, Chakravarty K, et al. BSR and BHPH guideline for the management of adults with ANCA-associated vasculitis. *Rheumatology (Oxford).* 2014;53:2306. <http://rheumatology.oxfordjournals.org/content/53/12/2306/suppl/DC1> Published April 11, 2014. Accessed August 06, 2015.
- Kaldas A, Warraich I, Prabhakar SS. ANCA associated glomerulonephritis—an in-depth review. *J Nephrol Ther.* 2013;4:doi:10.4172/2161-0959.1000147.
- Wu YY, Hsu TC, Chen TY, et al. Proteinase 3 and dihydrolipoamide dehydrogenase (E3) are major autoantigens in hepatitis C virus (HCV) infection. *Clin Exp Immunol.* 2002;128:347-352.
- Kemna MJ, Damoiseaux J, Austen J, et al. ANCA as a predictor of relapse: useful in patients with renal involvement but not in patients with nonrenal disease. *J Am Soc Nephrol.* 2015;26:537-542.

\* The CPT codes provided are based on AMA guidelines and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.