


The Clinical Insights Companion, an eClinicalWorks® app

A “how-to” guide to launch your new app as part of your daily EHR workflows


Easy, intuitive steps to begin using your Clinical Insights Companion app today

1. Access the app via the eClinicalWorks® (eCW) patient dashboard


A) View the Patient Profile



B) Select “Apps—Quest Diagnostics” to launch the Clinical Insights Companion (CIC) app directly in your workflow




C) Click on “I agree” to accept the disclaimer/patient consent




If you click to cancel, the “Not verified” messaging appears

If you click on “I agree,” you will see the options as outlined below


2. Access multiple views of the patient’s diagnostic history based on your workflow preferences



View list of results by different time frames




Open the **Test groups** drop-down to filter results by clinical category




Expanding the **Test groups** drop-down shows results in your selected clinical category

2. Access multiple views of the patient diagnostic history based on your workflow preferences (continued)




Search for a specific test by test name or code




Filter to only out-of-range (high/low) results using the toggle view


3. View results using either the List view or the Chart view



View result trends and changes over time using the Chart view




View diagnoses, prescriptions, and allergies without leaving the interface



View encounter details, including diagnosis, provider, and date of service


4. View specific test details



Hover over any result to view detailed test information

You can hover over any test result in both List and Chart views, giving you quick access to key test information without leaving your current screen. Such details may include

- Reported date and time
- Reference range
- Ordering physician
- Specimen ID
- Data source



Print results directly from the detail view

Tips

- Use the Chart view to visualize results and trends over time
- Use the Up and Down arrows to highlight distinct data points to easily identify high, low, and normal result values
- Switch back to the List view for detailed result values



Quest is dedicated to improving the health of your patients and your organization. Visit us at [QuestDiagnostics.com](https://www.questdiagnostics.com) for further information on the Quest CIC app for eCW and how it consolidates patient diagnostic history into one easy-to-navigate, time-sequenced view within your EHR.