


The Clinical Insights Companion, an eClinicalWorks® app

A “how-to” guide to launch your new app as part of your daily EHR workflows


Easy, intuitive steps to begin using your Clinical Insights Companion™ (CIC) app today

1. Access the app via the eClinicalWorks® (eCW) patient dashboard


A) View the Patient Profile



B) Select “Apps—Quest Diagnostics” to launch the Clinical Insights Companion (CIC) app directly in your workflow




C) Click on “I agree” to accept the disclaimer/patient consent




If you click to cancel, the “Not verified” messaging appears

If you click on “I agree,” you will see the options as outlined below


2. Access multiple views of the patient lab result history with data from the patient's chart based on your workflow preferences



View list of results by different time frames



Open the Test groups drop-down to filter results by clinical category




Expanding the Test groups drop-down shows results in your selected clinical category

2. Access multiple views of the patient lab result history with data from the patient's chart based on your workflow preferences (continued)

Search for a specific test by test name or code

Filter to only **out-of-range** (high/low) results using the toggle view

3. View results using either the List view or the Chart view




View result trends and changes over time using the **Chart** view

View **diagnoses, prescriptions, and allergies** without leaving the interface

View **encounter details**, including diagnosis, provider, and date of service

4. View specific test details



Hover over any result to view detailed test information

You can hover over any test result in both List and Chart views, giving you quick access to key test information without leaving your current screen. Such details may include

- Reported date and time
- Reference range
- Ordering physician
- Specimen ID
- Data source

Print results directly from the detail view

Tips

- Use the Chart view to visualize results and trends over time
- Use the Up and Down arrows to highlight distinct data points to easily identify high, low, and normal result values
- Switch back to the List view for detailed result values



Quest is dedicated to improving the health of your patients and your organization. Visit us at **QuestDiagnostics.com** for further information on the Quest CIC app for eCW and how it consolidates patient lab result history with data from the patient's chart into one easy-to-navigate, time-sequenced view within your EHR.

QuestDiagnostics.com

Quest®, Quest Diagnostics®, any associated logos, and all associated Quest Diagnostics registered or unregistered trademarks are the property of Quest Diagnostics. All third-party marks—® and ™—are the property of their respective owners. © 2026 Quest Diagnostics Incorporated. All rights reserved. FAQ14301 2/2026

