

5 Panel DOT look-alike

Unit code 65105N - Substance abuse panel

Drug class	Initial test level	Confirmatory test level
Amphetamines Amphetamine Methamphetamine	500 ng/mL	250 ng/mL 250 ng/mL
Cocaine metabolites	150 ng/mL	100 ng/mL
Marijuana metabolites	50 ng/mL	15 ng/mL
MDA-analogues MDA MDMA	500 ng/mL	250 ng/mL 250 ng/mL
Opiates Codeine Morphine	2000 ng/mL	2000 ng/mL 2000 ng/mL
6-Acetylmorphine	10 ng/mL	10 ng/mL
Opiates (semi-synthetic) Hydromorphone Hydrocodone	300 ng/mL	100 ng/mL 100 ng/mL
Oxycodones Oxymorphone Oxycodone	100 ng/mL	100 ng/mL 100 ng/mL
Phencyclidine	25 ng/mL	25 ng/mL
TestSure™		



10 Panel DOT look-alike

Unit code 65190N - Substance abuse panel

Drug class	Initial test level	Confirmatory test level
Amphetamines Amphetamine Methamphetamine	500 ng/mL	250 ng/mL 250 ng/mL
Barbiturates Amobarbital Butalbital Pentobarbital Phenobarbital Secobarbital	300 ng/mL	200 ng/mL 200 ng/mL 200 ng/mL 200 ng/mL 200 ng/mL
Benzodiazepines Alprazolam metabolite Oxazepam	300 ng/mL	200 ng/mL 200 ng/mL
Cocaine metabolites	150 ng/mL	100 ng/mL
Marijuana metabolites	50 ng/mL	15 ng/mL
Methadone	300 ng/mL	200 ng/mL
Methaqualone	300 ng/mL	200 ng/mL
MDA-analogues MDA MDMA	500 ng/mL	250 ng/mL 250 ng/mL
Opiates Codeine Morphine	2000 ng/mL	2000 ng/mL 2000 ng/mL
6-Acetylmorphine	10 ng/mL	10 ng/mL
Opiates (semi-synthetic) Hydromorphone Hydrocodone	300 ng/mL	100 ng/mL 100 ng/mL
Oxycodones Oxymorphone Oxycodone	100 ng/mL	100 ng/mL 100 ng/mL
Phencyclidine	25 ng/mL	25 ng/mL
Propoxyphene	300 ng/mL	200 ng/mL
TestSure™		

Employer Solutions | 1.800.877.7484

Commonly ordered urine unit codes



5 Panel pre 2010 DOT look-alike

Unit code 35105N - Substance abuse panel

Drug class	Initial test level	Confirmatory test level
Amphetamines Amphetamine Methamphetamine	1000 ng/mL	500 ng/mL 500 ng/mL
Cocaine metabolites	300 ng/mL	150 ng/mL
Marijuana metabolites	50 ng/mL	15 ng/mL
Opiates Codeine Morphine	2000 ng/mL	2000 ng/mL 2000 ng/mL
Phencyclidine	25 ng/mL	25 ng/mL
TestSure™		



10 Panel pre 2010 DOT look-alike

Unit code 35190N - Substance abuse panel

Drug class	Initial test level	Confirmatory test level
Amphetamines Amphetamine Methamphetamine	1000 ng/mL	500 ng/mL 500 ng/mL
Barbiturates Amobarbital Butalbital Pentobarbital Phenobarbital Secobarbital	300 ng/mL	200 ng/mL 200 ng/mL 200 ng/mL 200 ng/mL 200 ng/mL
Benzodiazepines Alprazolam metabolite Oxazepam	300 ng/mL	200 ng/mL 200 ng/mL
Cocaine metabolites	300 ng/mL	150 ng/mL
Marijuana metabolites	50 ng/mL	15 ng/mL
Methadone	300 ng/mL	200 ng/mL
Methaqualone	300 ng/mL	200 ng/mL
Opiates Codeine Morphine	2000 ng/mL	2000 ng/mL 2000 ng/mL
Phencyclidine	25 ng/mL	25 ng/mL
Propoxyphene	300 ng/mL	200 ng/mL
TestSure™		



Expanded opiate & benzodiazepine

Unit code 32093N - Substance abuse panel

Drug class	Initial test level	Confirmatory test level
Amphetamines Amphetamine Methamphetamine	1000 ng/mL	500 ng/mL 500 ng/mL
Barbiturates Amobarbital Butalbital Pentobarbital Phenobarbital Secobarbital	300 ng/mL	200 ng/mL 200 ng/mL 200 ng/mL 200 ng/mL 200 ng/mL
Benzodiazepines Alprazolam metabolite Oxazepam Flurazepam metabolite Lorzapam Nordiazepam Temazepam Triazolam metabolite	300 ng/mL	200 ng/mL 200 ng/mL 200 ng/mL 200 ng/mL 200 ng/mL 200 ng/mL 200 ng/mL
Cocaine metabolites	300 ng/mL	150 ng/mL
Marijuana metabolites	50 mg/mL	15 ng/mL
Methadone	300 ng/mL	200 ng/mL
MDA-analogues MDA MDMA MDEA	500 ng/mL	250 ng/mL 250 ng/mL 250 ng/mL
Opiates Codeine Morphine Hydromorphone Hydrocodone	300 ng/mL	300 ng/mL 300 ng/mL 300 ng/mL 300 ng/mL
6-Acetylmorphine	10 ng/mL	10 ng/mL
Oxycodones Oxymorphone Oxycodone	100 ng/mL	100 ng/mL 100 ng/mL
Phencyclidine	25 ng/mL	25 ng/mL
Propoxyphene	300 ng/mL	200 ng/mL
TestSure™		

Employer Solutions | **1.800.877.7484**

Commonly ordered urine unit codes