


Express Results™ Online

An instant, urine drug testing solution designed to deliver fast, dependable results

Employers depend on reliable drug testing results to help them make informed hiring and staffing decisions. When turnaround time is of the essence, Express Results™ Online provides a fast, accurate, electronically-driven drug testing solution. Negative drug screen results are available within minutes from the time of collection for up to 11 different drugs. Employers can order the drug test, locate a collection site, email the donor an authorization form, access real-time status updates throughout the process and get the final results — all simply and conveniently online.

- Fast results available in minutes
- Detects a wide range of drugs
- Collection sites staffed with trained collectors
- Improved efficiency using technology
- Secure systems ensure confidentiality and privacy

A fast, accurate process with results available in minutes



Screen for up to 11 drugs

Express Results *Online* screens for up to 11 different drugs including barbiturates, benzodiazepines, amphetamine, cocaine, marijuana, methadone, methamphetamine, opiates, oxycodone, PCP and propoxyphene. Employers can choose a pre-selected combination of drugs or configure a custom panel.

Specimen validity testing is performed on every Express Results *Online* specimen in order to detect an invalid, adulterated or substituted sample. We check creatinine, oxidants (e.g. nitrites) and pH levels in addition to checking the temperature of every urine sample at the time of collection.

| Drug | Cut-off |
|---------------------|-------------|
| Amphetamine | 1,000 ng/mL |
| Methamphetamine | 1,000 ng/mL |
| Cocaine | 300 ng/mL |
| Opiates | 2,000 ng/mL |
| Phencyclidine (PCP) | 25 ng/mL |
| Barbiturates | 300 ng/mL |
| Benzodiazepines | 300 ng/mL |
| Oxycodone | 100 ng/mL |
| Marijuana | 50 ng/mL |
| Methadone | 300 ng/mL |
| Propoxyphene | 300 ng/mL |

Get accurate results within minutes

Express Results *Online* empowers employers to easily manage the drug testing process electronically with remarkable speed. Negative results are typically reported within minutes of the drug test collection. Non-negative specimens are transported to one of our four, SAMHSA-certified laboratories for confirmation testing. The laboratory typically reports the final result within 24 to 48 hours of the specimen's arrival.

Simplify the collection experience

We provide unparalleled access to a combined nationwide network of company-owned and operated Quest Diagnostics Patient Service Centers and preferred, third-party collection sites for drug testing. You can locate the nearest collection site, set up an appointment and print maps — all can be done from the convenience of your computer. Every collector is trained to perform an Express Results *Online* test and follows a step-by-step, online protocol that helps to ensure correct, complete information is captured. Because a certified collector reads the drug test consistently in a controlled environment, the integrity of the process is maintained, giving employers added peace of mind and confidence in the final, reported result.

To learn more about Express Results *Online* and how it can benefit your organization, contact your Quest Diagnostics representative today.

We're there when you need us™
EmployerSolutions.com | 800.877.7484